



Biosolve Chimie
20 Rue Roger Husson, 57260 Dieuze, France.

Tel: +33 3 878 675 80 Fax: +33 3 878 675 89
Email: info@biosolve-chimie.com

Certified ISO 9001:2015

Biosolve BV
Kerkhofstraat 21, 5554HG Valkenswaard, the Netherlands.

Tel: +31-(0)40-2071300 Fax: +31-(0)40-204 85 37
Email: info@biosolve-chemicals.com

Internet website: www.biosolve-chemicals.com

CERTIFICATE OF COMPLIANCE/ANALYSIS

Product : Dimethylsulfoxide Headspace 1L

Synonym: DMSO, Methyl sulfoxide.



Article no : 0004477501BS

CAS 67-68-5

M.P. 18.4 °C

Lot#: 0010009759

EC NO 200-664-3

B.P. 189 °C

Grade: Headspace

FW C2H6OS

MW 78.13

DENSITY 1.1

Test	Specification	Result
Purity (GC)	>= 99.99 %	99.99 %
Identity (IR)	Complying	Complying
Water (KF)	<= 0.05 %w/w	< 0.01 %w/w
Color (APHA Pt-Co)	<= 10	5
GC-HS (Each class 1 solvent)	<= 1 µg/g	< 1 µg/g
GC-HS (Each class 2 solvent)	<= 10 µg/g	< 10 µg/g
GC-HS (Each class 3 solvent)	<= 50 µg/g	< 50 µg/g
Headspace test for O.V.I's	Complying	Complying

Suitable as solvent for the analysis of residual solvents of ICH classes 1,2 and 3 according to Ph Eur and USP.
Formic acid, acetic acid, dimethylacetamide and cumene are not detected.

Manufacturing Date : 12/2022

Expiration Date : 12/2025

TERMS OF SALE:

SELLERS ONLY OBLIGATION SHALL BE TO REPLACE SUCH QUANTITY OF THE PRODUCT PROVEN TO BE DEFECTIVE. SELLER SHALL NOT BE LIABLE FOR ANY INJURY, LOSS OR DAMAGE, DIRECT OR CONSEQUENTIAL, ARISING OUT OF THE USE, MISUSE, OR THE INABILITY TO USE THE PRODUCT. USER SHALL DETERMINE THE SUITABILITY OF THE PRODUCT FOR ITS INTENDED USE AND SHALL ASSUME ALL LIABILITY IN CONNECTION THEREWITH.

Store at room temperature in a shaded and ventilated place.

Our COA's are computer printouts and therefore valid without signature.

SDS is available on our website.

Y.Shalom Ph.D.
QC Manager