



**Biosolve Chimie**  
20 Rue Roger Husson, 57260 Dieuze, France.

Tel: +33 3 878 675 80 Fax: +33 3 878 675 89  
Email: info@biosolve-chimie.com

**Certified ISO 9001:2015**

**Biosolve BV**  
Kerkhofstraat 21, 5554HG Valkenswaard, the Netherlands.

Tel: +31-(0)40-2071300 Fax: +31-(0)40-204 85 37  
Email: info@biosolve-chemicals.com

Internet website: www.biosolve-chemicals.com

## CERTIFICATE OF COMPLIANCE/ANALYSIS

Product : Methanol absolute HPLC 4 L

Synonym: Methyl alcohol.



Article no : 0013680604BS

CAS 67-56-1

M.P. -98 °C

Lot#: 0010010851

EC NO 200-659-6

B.P. 64.7 °C

Grade: HPLC

FW CH4O

n20/D1.329

MW 32.04

DENSITY 0.791

UN 1230,3,II

Test	Specification	Result
Clear colorless liquid	Complying	Complying
Acidity (as Acetic acid)	<= 0.002 %	0.001 %
Color (APHA Pt-Co)	<= 10	5
Assay (GC on anhydrous basis)	>= 99.9 %	100.0 %
Grad. elution H.Peak at 254nm	<= 0.002 AU	0.001 AU
Residue after evaporation	<= 0.0003 %w/w	< 0.0001 %w/w
Water (KF)	<= 0.05 %w/w	0.02 %w/w
UV Trans. 210nm	>= 30 %	34 %
UV Trans. 235nm	>= 85 %	86 %
UV Trans. 260nm	>= 98 %	99 %

Filtered through 0.2um, filled under inert gas.

Manufacturing Date : 02/2023

Expiration Date : 02/2028

#### TERMS OF SALE:

SELLERS ONLY OBLIGATION SHALL BE TO REPLACE SUCH QUANTITY OF THE PRODUCT PROVEN TO BE DEFECTIVE. SELLER SHALL NOT BE LIABLE FOR ANY INJURY, LOSS OR DAMAGE, DIRECT OR CONSEQUENTIAL, ARISING OUT OF THE USE, MISUSE, OR THE INABILITY TO USE THE PRODUCT. USER SHALL DETERMINE THE SUITABILITY OF THE PRODUCT FOR ITS INTENDED USE AND SHALL ASSUME ALL LIABILITY IN CONNECTION THEREWITH.

Store at room temperature in a shaded and ventilated place.

Our COA's are computer printouts and therefore valid without signature.

SDS is available on our website.

Y.Shalom  
Ph.D. QC