



TOMO® IHC Adhesive Glass Slides

- Intended Use -

Matsunami TOMO® IHC Adhesive Glass Microscope Slides are recommended for use on Ventana® Benchmark staining platform and automated IHC platforms.

- General Information -

Matsunami TOMO® hydrophilic adhesion slides offer superior tissue adhesion, staining performance and slide surface fluidics. Proprietary slide coating ensures consistent performance even after exposure to heat stress and humid environmental conditions. Optimized for use with thermal printers. Matsunami microscope slides are manufactured in Japan using optically clear, exceptionally flat, low fluorescence glass and finished with 90° and clipped corners. ISO 13485 Certified. IVD CE marked.

- Specifications -

Case Quantity	1000	<i>Slides</i>
Case Width	6.625	<i>Inches</i>
Case Length	9.75	<i>Inches</i>
Case Height	4.125	<i>Inches</i>
Case Weight	11.55	<i>Lbs</i>
Slide Width	25	<i>Milimeters</i>
Slide Height	75	<i>Milimeters</i>
Slide Thickness	1.0	<i>Milimeters</i>
Paint Height	20.14	<i>Milimeters (±0.15mm)</i>
Wettability	Hydrophilic Surface	
Corner Options	45° or 90°	
Color Options	Blue, Green, Pink, White, Yellow	
Origin Country	Japan	

- Packaging -

Slides are packaged 100 per box, 10 boxes per case (1000 slides)

Catalog#	Color	Corner
TOM-1190	White	90°
TOM-11	White	45°
TOM-12	Yellow	45°
TOM-13	Pink	45°
TOM-14	Blue	45°
TOM-15	Green	45°

